

# Particle Tracker LabView Code

Wednesday, July 01, 2015 1:13 PM

- . **Step 1:** Open the video file in ImageJ (\*.avi)
- . **Step 2:** Rotate the video such that the channel is horizontal.
- . **Step 3:** Crop the video using a rectangle 1 pixel wide and as high as the channel where the laser is positioned.
- . **Step 4:** Save the cropped video as image sequence (text).
- . **Step 5:** Use the LabView code to analyze the images.
  - Set the interval at 20-25.
  - Set the Intensity threshold base on the normalized intensity plot.
  - Pixel/um is ~0.74 for 20x objective.


## Code Versions:

V1	<a href="file:///C:/Users/mehranh/Documents/Results/Inertial Focusing/Cell Intensity Measurement.exe">file:///C:/Users/mehranh/Documents/Results/Inertial Focusing/Cell Intensity Measurement.exe</a>	First version including threshold and histogram.
V1.1	<a href="file:///C:/Users/mehranh/Documents/Results/Inertial Focusing/Cell Intensity Measurement 1.1/Cell Intensity Measurement.exe">file:///C:/Users/mehranh/Documents/Results/Inertial Focusing/Cell Intensity Measurement 1.1/Cell Intensity Measurement.exe</a>	Loop version of V1 plus pixel/um added
V1.2	<a href="file:///C:/Users/mehranh/Documents/Results/Inertial Focusing/Cell Intensity Measurement 1.2/Cell Intensity Measurement.exe">file:///C:/Users/mehranh/Documents/Results/Inertial Focusing/Cell Intensity Measurement 1.2/Cell Intensity Measurement.exe</a>	Added gating capabilities; by setting the gate center and half width, the percentage of events in the gate is shown.
V1.3	<a href="file:///C:/Users/mehranh/Documents/Results/Inertial Focusing/Cell Intensity Measurement 1.3/Cell Intensity Measurement.exe">file:///C:/Users/mehranh/Documents/Results/Inertial Focusing/Cell Intensity Measurement 1.3/Cell Intensity Measurement.exe</a>	Gating parameters changed to gate center and width (used to be half width)
V1.4	<a href="file:///C:/Users/mehranh/Documents/Results/Inertial Focusing/Cell Intensity">file:///C:/Users/mehranh/Documents/Results/Inertial Focusing/Cell Intensity</a>	Gate bounds are now shown on the

V1.4	< <a href="file://C:\Users\mehranh\Documents\Results\Inertial Focusing\Cell Intensity Measurement 1.4\Cell Intensity Measurement.exe">file://C:\Users\mehranh\Documents\Results\Inertial Focusing\Cell Intensity Measurement 1.4\Cell Intensity Measurement.exe</a> >	Gate bounds are now shown on the histogram as dotted red lined
V1.5	< <a href="file://C:\Users\mehranh\Documents\Results\Inertial Focusing\Cell Intensity Measurement 1.5\Cell Intensity Measurement.exe">file://C:\Users\mehranh\Documents\Results\Inertial Focusing\Cell Intensity Measurement 1.5\Cell Intensity Measurement.exe</a> >	Added plot for intensity profile of each frame.
V1.5.1	< <a href="file://C:\Users\mehranh\Documents\Results\Inertial Focusing\Cell Intensity Measurement 1.5\Cell Intensity Measurement.exe">file://C:\Users\mehranh\Documents\Results\Inertial Focusing\Cell Intensity Measurement 1.5\Cell Intensity Measurement.exe</a> >	Frame picker line moves to the center of the new range
V1.6	< <a href="file://C:\Users\mehranh\Documents\Results\Inertial Focusing\Cell Intensity Measurement 1.5\Cell Intensity Measurement.exe">file://C:\Users\mehranh\Documents\Results\Inertial Focusing\Cell Intensity Measurement 1.5\Cell Intensity Measurement.exe</a> >	Added plot for intensity profile of each frame plus a Gaussian fit and peaks found by peakfind.
V1.7	 Cell Intensity ...	<ul style="list-style-type: none"> <li>• Removed white space at end of intensity data</li> <li>• Intensity data is normalized</li> <li>• Cell peaks are displayed on plot</li> <li>• Added histogram for cell arrival time delta</li> </ul>
V1.7.1	 Cell Intensity ...	<ul style="list-style-type: none"> <li>• Added support of 2D array files</li> <li>• FWHM of Gaussian fit reported (for visualization only)</li> </ul>
V2.0	 Cell Intensity ...	<ul style="list-style-type: none"> <li>• New tabbed interface</li> <li>• New normalization routine. Average frames w/o peaks</li> <li>• Gaussian fit identifies cells</li> <li>• Cell width histogram</li> </ul>

## V1.6 sample frames

Wednesday, July 15, 2015 2:12 PM



Strategy for next version:

- Calculate the baseline by averaging the minimum intensity, and the two ends of each file.
- For frames with "two peaks detected by peakfind" and "maximum value = 255" like example #6, a Gaussian fit is performed by ignoring the flat part of the profile. Then, the FWHM and location of the peak is extracted.
- For frames with "maximum value in less than 255" or "only one max value = 255" like example #4 and #6, the FWHM is calculated by marching to the sides of the peak down to the half magnitude. The location is also determined based on the midpoint of the Half Magnitude line.

Intensity fitting strategy

Find Baseline by  
averaging Min and  
first and last pixel  
intensities

